

National Health and Nutrition Examination Survey

2017-March 2020 Data Documentation, Codebook, and Frequencies

Urine Pregnancy Test (P_UCPREG)

Data File: P_UCPREG.xpt

First Published: June 2021

Last Revised: NA

Component Description

The NHANES program suspended field operations in March 2020 due to the coronavirus disease 2019 (COVID-19) pandemic. As a result, data collection for the NHANES 2019-2020 cycle was not completed and the collected data are not nationally representative. Therefore, data collected from 2019 to March 2020 were combined with data from the NHANES 2017-2018 cycle to form a nationally representative sample of NHANES 2017-March 2020 pre-pandemic data. These data are available to the public. Please refer to the Analytic Notes section for more details on the use of the data.

A urine pregnancy test was performed on menstruating female survey participants 8 years and older. All positive test results were used to exclude pregnant women from the Dual Energy X-ray Absorptiometry component at the mobile examination center (MEC).

Eligible Sample

Examined female participants aged 12–59 years, and menstruating females aged 8–11 years, in the NHANES 2017-March 2020 pre-pandemic sample, were eligible. However, due to disclosure risks, only females 20-44 years of age have urine pregnancy results in this file.

Description of Laboratory Methodology

The Icon 25 hCG test kit (Beckman Coulter) is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine or serum to aid in the early detection of pregnancy. The test utilizes a combination of monoclonal and polyclonal antibodies to selectively detect elevated levels of hCG in urine or serum.

Laboratory Quality Assurance and Monitoring

Urine specimens were processed in the NHANES mobile examination centers.

Detailed instructions on specimen collection and processing are discussed in the [2017-2018](#) and [2019-2020 NHANES Laboratory Procedures Manuals](#) (LPMs).

The NHANES quality assurance and quality control (QA/QC) protocols meet the 1988 Clinical Laboratory Improvement Act mandates. Detailed QA/QC instructions are discussed in the NHANES LPMs.

Mobile Examination Centers (MECs)

Laboratory team performance is monitored using several techniques. NCHS and contract consultants use a structured competency assessment evaluation during visits to evaluate both the quality of the laboratory work and the QC procedures. Each laboratory staff member is observed for equipment operation, specimen collection and preparation; testing procedures and constructive feedback are given to each staff member. Formal retraining sessions are conducted annually to ensure that required skill levels were maintained.

Data Processing and Editing

The data were reviewed. Incomplete data or improbable values were sent to the performing laboratory for

confirmation.

Analytic Notes

The COVID-19 pandemic required suspension of NHANES 2019-2020 field operations in March 2020 after data were collected in 18 of the 30 survey locations in the 2019-2020 sample. Data collection was cancelled for the remaining 12 locations. Because the collected data from 18 locations were not nationally representative, these data were combined with data from the previous cycle (2017-2018) to create a 2017-March 2020 pre-pandemic data file. A special weighting process was applied to the 2017-March 2020 pre-pandemic data file. The resulting sample weights in the present file should be used to calculate estimates from the combined cycles. These sample weights are not appropriate for independent analyses of the 2019-2020 data and will not yield nationally representative results for either the 2017-2018 data alone or the 2019-March 2020 data alone. Please refer to the NHANES website for additional information for the NHANES 2017-March 2020 pre-pandemic data, and for the previous 2017-2018 public use data file with specific weights for that 2-year cycle.

Refer to the [2017-2018](#) and [2019 - 2020 Laboratory Data Overview](#) documents for general information on NHANES laboratory data.

There are over 800 laboratory tests performed on NHANES participants. However, not all participants provided biospecimens or enough volume for all the tests to be performed. The specimen availability can also vary by age or other population characteristics. Analysts should evaluate the extent of missing data in the dataset related to the outcome of interest as well as any predictor variables used in the analyses to determine whether additional re-weighting for item non-response is necessary.

Please refer to the NHANES [Analytic Guidelines](#) and the on-line NHANES [Tutorial](#) for further details on the use of sample weights and other analytic issues.

Demographic and Other Related Variables

The analysis of NHANES laboratory data must be conducted using the appropriate survey design and demographic variables. The [NHANES2017-March 2020 Pre-pandemic Demographics File](#) contains demographic data, health indicators, and other related information collected during household interviews as well as the sample design variables. The recommended procedure for variance estimation requires use of stratum and PSU variables (SDMVSTRA and SDMVPSU, respectively) in the demographic data file.

This laboratory data file can be linked to the other NHANES data files using the unique survey participant identifier (i.e., SEQN).

Detection Limits

Since this data is reported as qualitative data, the use of lower LLOD is not applicable.

Codebook and Frequencies

SEQN - Respondent sequence number

Variable Name:	SEQN
SAS Label:	Respondent sequence number
English Text:	Respondent sequence number.
Target:	Females only 20 YEARS - 44 YEARS

URXPREG - Urine Pregnancy Result

Variable Name: URXPREG
SAS Label: Urine Pregnancy Result
English Text: Urine Pregnancy Result
Target: Females only 20 YEARS - 44 YEARS

Code or Value	Value Description	Count	Cumulative	Skip to Item
1	Positive	84	84	
2	Negative	1630	1714	
3	Not Done	33	1747	
4	Invalid	0	1747	
.	Missing	3	1750	